

# **CAMELID VHH** FRAGMENT



**Expert in Antibody Engineering** 

From Research to Discovery



ISO 9001 : 2015 certified

## **PROTOCOL OVERVIEW**



### **ABOUT VHH**

Camelids produce conventional antibodies made of 2 heavy chains and 2 light chains bound together with disulphide bonds in a Y shape. However, they also produce a unique subclasse of immunoglobulin G, also known as heavy chain antibody (HC-Ab).

These antibodies are made of only two heavy chains without the CH1 region, and an antigen binding domain at their N-terminus called VHH.

The unique feature of hclgG is the capacity of their monomeric antigen binding regions (VHH) to bind antigens with specificity, affinity and especially diversity that are comparable to conventional antibodies without the need of pairing with another region. Using molecular engeenering VHH fragments can be generated.

### COVALAB, YOUR PARTNER IN ANTIBODY ENGINEERING SINCE 1995

### MADE IN FRANCE - ISO 9001 : 2015 CERTIFIED



### WHY USE VHH FRAGMENTS RATHER THAN CONVENTIONAL ANTIBODIES?

Small & accessible From 12 to 15 kDa - Bind the hidden epitope - Reach hardly accessible antigens

**Penetration** Quicker - Target intracellular proteins - Precise staining of tissue sections

> Soluble & stable Thermostable up to 70°C / 158°F - Resist to extreme pH

Easy to manufacture Easy to manufacture in yeast or microbial systems

> Specificity & affinity Nano to picomolar affinity

Clearance Rapid blood elimination - Able to pass renal filter

### **DETAILED PROTOCOLE**



### **IMMUNOGEN SELECTION**



Protein



**Microorganism** Inactivated virus Inactivated yeast Inactivated bacteria



Peptide PTM



Hapten Carbohydrates Chemicals & Toxins Natural or modified nucleotides





Other Cell extracts Tissus extracts Others (please contact us)



### LLAMA IMMUNIZATION

According to our exclusive protocol, the immunization last from 53 to 88 days, according to the results, with up to 4 injections and 2 test sera. The immunoreactivity is monitored by the detection method developed for the project.



All the experiments are carried out by experienced and authorised personnel according to H&S procedures, established in accordance to the French legislation.

### **YOU RECEIVE**

according to the immunoreactivity.

Serua to run tests in your specific conditions and adapt the protocol

### **GUARANTEE**

- ß
- The immunized animal is kept in the animal facility to perform additional total RNA extraction if necessary.



### PHAGE LIBRARY CONSTRUCTION



**PBMCs isolation and RNA extraction** 

Peripheral blood mononuclear cells (PBMCs) are isolated from the final serum and total RNA is extracted.

### cDNA synthesis and VHH amplification

RT-PCR and VHH antibody fragment amplification Specific primers are design to perform retro-transcription and amplification of VHH antibody fragments.

### Phage vector cloning

VHH coding genes are integrated into a phage vector which is then used to transform competent bacteria for phage library production.



### PANNING

Several selection rounds are required for target specific recombinant phage enrichment.





### **HIGH SCALE PRODUCTION**

The cDNAs coding for VHH of interest can then be used for high-scale production. The VHH gene are subcloned into our expression vector. VHH are then produced after transformation and growth of competent bacteria followed by affinity chromatography purification.



PCR amplification of the selected

VHH gene.

100000 C





→ V<sub>µ</sub>H

Inserts are then integrated into a high yield expression vector which is used to transform competent bacteria. Specific VHH are purified from cells lysis supernatants by affinity chromatography.

### YOU RECEIVE

 VHH preparations purified from bacteria cell culture medium by antigen-affinity chromatography

### **GUARANTEE**

 All cDNAs coding for specific VHH are kept frozen at -80°C in our laboratory



# SCIENTIFIC SUPPORT EXPERTISE HIGH QUALITY REACTIVITY FIFXIBI SINCE 19



### Covalab S.A.S.

1B rue Jacques Monod 69500 Bron France

► +33 (0)4 37 65 42 30
☑ enquiries@covalab.com

www.covalab.com



ISO 9001:2015 BUREAU VERITAS Certification

