

Antibody Development Services

Recombinant Nanobodies

www.covalab.com

WHY USE NANOBODIES?



Immunoglobulins are very powerful tools which are used in many lab experiments. But over the last few years, their size has progressively become an issue for some applications. Many solutions have been engineered to overcome this problem, such as Fab or single chain (ScFv) fragments.



Camelids \lg_2 and \lg_3 possess only 2 heavy chains composed of two constant regions and a very versatile antigen binding domain called VHH. Using molecular engineering, high affinity VHH fragments of only 15 kDa called Nanobodies, which are more stable and soluble than conventional \lg_6 , can now be produced in high scale, thus allowing to reach hardly accessible antigens and penetrate tissues easier.

WHICH TYPE OF IMMUNOGEN?

PROTEIN

To perform all required steps of the immunization protocol, which includes immunization, development of the screening method and immunoreactivity monitoring, we require a few milligrams of the immunogen (recommended concentration: 2 mg/ml), which can be as:



Protein

In solution Lyophilized Within polyacrylamide gel



Microorganism

Inactivated virus Inactivated bacteria Inactivated yeast



Other

Cell or tissue extracts Others (please contact us)

PEPTIDE

POST-TRANSLATIONAL MODIFICATION

Design

We use our in-house expertise to design the most relevant peptides to ensure the specificity of the nanobody.

We analyze **physicochemical properties** to predict immunogenicity and tertiary structure among others.

We ensure **minimal cross-reactivity** with non-relevant proteins by analyzing their homology with other sequences using alignment tools.

Synthesis

Peptides are synthesized using solid phase Fmoc chemistry, and analyzed by HPLC and mass spectrometry.

Length: 15 residues in average

Quantity: 20 to 25 mg

When developing an **anti-PTM** antibody, we synthesize one modified peptide as well as **one unmodified control peptide** to ensure the specificity of the resulting antibody for the modification.

HAPTENS

Conjugation

Conjugation is usually carried out on a **cysteine residue** (which can be added to the *N- or C-terminus of the peptide if required*), but other techniques can be used if necessary. After synthesis, peptides are conjugated with one of the following carrier proteins:

KLH Keyhole limpet hemocyanin



BSA Bovine serum albumin



Carbohydrates



Chemicals & Toxins



Natural or modified nucleotides

DETAILED PROTOCOL

IMMUNIZATION

The first step of the procedure consists in injecting the immunogen into the host llama to generate an immune response, according to our 88-day exclusive protocol dedicated to this animal. Thanks to our 20-year experience and know-how, the number and frequency of injections are optimized real-time to ensure the best immunoreactivity / quantity of immunogen ratio.



Each test bleed is assayed according to the detection method developed for the project. Our immunization protocols are adjustable according to the immunoreactivity of the sera. In case of weak immune response, additional injections and/or bleeds can be requested.

All experiments are undertaken by experienced and authorized staff following Health and Safety procedures, established according to the French legislation governing the use of animals in experiments.

You receive:



 Test bleeds so that you can run tests in your specific conditions in order to adapt the immunization protocol if the immunoreactivity is not as high as expected. • The immunized animal is kept in the animal facility to perform additional total RNA extraction if necessary.

PHAGE LIBRARY CONSTRUCTION

After confirmation of the secretion of antibodies into the bloodstream according to both our ELISA tests and your results, we proceed to the extraction of total RNA from PBMCs. Using suitable primers, retro-transcription is performed to obtain a library of cDNAs which are then inserted into phage vectors. Transformation into bacteria allows the expression of corresponding proteins at the surface of the phages.



Peripheral blood mononuclear cells (PBMCs) are isolated from the final bleed and total RNA is extracted. Specific primers are designed to perform retro-transcription within the region of the genome coding for VHH and second nested PCR allows cDNA amplification. • cDNAs are integrated into our phage vector which is then used to transform competent bacteria for phage production. Phages are produced by the bacteria and released into the culture medium upon lysis.

Our guarantee:

PANNING



Several rounds of selection are necessary to isolate the phages which carry the cDNA coding for the VHH of interest. By immobilizing the antigen in microtiter plate wells, only phages which can bind the antigen are retained while the others are removed. Transduction of bacteria using those phages allows enrichment and amplification of specific VHH-coding cDNAs.



HIGH SCALE PRODUCTION



The cDNAs coding for VHH of interest can then be used for high-scale production. The nucleic acid sequences are extracted from phage genome and are subcloned into our expression vector. Nanobodies are then produced after transformation and growth of competent bacteria followed by affinity chromatography purification.



PROTOCOL OVERVIEW





Services

From project design to antibody production

Custom projects

Protocols adapted to your requirements

Technical support

Advice by experts in immunology & molecular biology



Approved Health & Safety procedures



www.covalab.com

Covalab S.A.S. — 11 avenue Albert Einstein 69100 Villeurbanne — France Phone: +33 (0) 437 654 230 — Fax: +33 (0) 437 289 416 E-mail: enquiries@covalab.com