

PRODUCT DATA SHEET

lab0162 Anti-Human GAM Immunoglobulins, Biotin conjugate antibody



Product specifications

Immunogen	IgGAM		
Form	Purified antibodies	Target	human
Raised in	Goat	Labelling	Biotin
Conjugate properties	The avidin-biotin system is a signal amplification system that allows the detection of quantities of antigens that are too small to be detected by direct enzymatic or fluorescent systems. Biotin labelled antibodies allow label changes and the use of multiple labels. Biotin-labelled antibodies have a long storage life. Moreover, by its small size, the biotin has a negligible effect on the biological activity of proteins. The presence of endogenous biotin could however, increase background.		
Applications	Antibodies are particularly suitable for use in the most sensitive procedures (ELISA, immuno-histochemistry, immuno-cytology and western blot). Their quality is guaranteed because a very stringent purification protocol is followed where non-specific and specific purifications alternate (affinity chromatography). Antibodies undergo a routine inspection of the cross reactions between target species. If necessary, the total elimination of cross-reactions can require an adsorption by affinity chromatography against the antigen whose recognition is not desired.		
Working dilutions	Optimal dilution should be determined by the end user The following are guidelines only : - ELISA : 1/1 000 to 1/10 000 - Immunoblot : 1/200 to 1/2 000 - Immunocytochemistry : 1/50 to 1/200 - Immunohistochemistry : 1/50 to 1/200		

Packaging specifications

Packaging	1 ml
Appearance	Liquid
Constituents	PBS, 50% Glycerol
Preservative	None
Storage	Short term storage : +4°C Long term storage : -20°C
Recommendations	Product in glycerol, do not freeze

Product type
Secondary antibody

Clonality
Polyclonal Antibody

Related product
opr0011 ; opr0002 ;
opr0008 ; opr0009

Warning
This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 04/10/06

version B

Covalab France
11, Avenue Albert Einstein
69100 Villeurbanne - France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab - UK Ltd
St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS - UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk