

PRODUCT DATA SHEET

mab0029 Monoclonal Antibody to Ochratoxin-A



Product Number/Size

mab0029-0 Lyophilized from 100 µl of ascitic fluid
mab0029-1 Lyophilized from 200 µl of ascitic fluid

Product specifications

Immunogen BSA-Ochratoxin-A

Form Ascitic fluid

Raised in Mouse

Clone 3C5

Isotype IgG1 / Lambda

Specificity Monoclonal antibody reacts specifically with Ochratoxin-A

Cross reactions Does not cross react with Ochratoxin-B

Working dilutions Optimal dilutions should be determined by the end user
The following are guidelines only :
- ELISA : 1/2 000 to 1/10 000

Packaging specifications

Appearance Lyophilized powder (glass vial)
Vial must be reconstituted as follow :
- mab0029-0 in 100 µl of distilled water
- mab0029-1 in 200 µl of distilled water

Preservative None

Storage Short term storage : +4°C
Long term storage : -20°C

References

Product type

Primary antibody

Clonality

Monoclonal Antibody

Related product

lab0252 ; lab0242 ;
lab0235 ; lab0246 ;
lab0257 ; lab0239 ;
lab0250 ...; opr0002;
opr0008 ; opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 05/02/2008

version C

Covalab France
11, Avenue Albert Einstein
69100 Villeurbanne - France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab - UK Ltd
St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS - UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk