

## PRODUCT DATA SHEET

### mab0050 Monoclonal Antibody to SNARE TI-VAMP Protein



#### Product Number/Size

mab0050-0 Lyophilized from 0,5 ml of tissue culture supernatant  
mab0050-1 Lyophilized from 1 ml of tissue culture supernatant

#### Product specifications

**Immunogen** GST-TI-VAMP recombinant protein

**Form** Tissue culture supernatant

**Raised in** Mouse

**Clone** 158.2

**Isotype** IgG1 / Kappa

**Specificity** Monoclonal antibody reacts specifically with the 25 kDa protein of human TI-VAMP protein (human, mouse, rat, canine)

**Cross reactions** Not yet tested

**Working dilutions** Optimal dilutions should be determined by the end user  
The following are guidelines only :  
- Immunohistology : 1/5 to 1/200  
- Western Blot : 1/50 to 1/500

#### Packaging specifications

**Appearance** Lyophilized powder (glass vial)  
Vial must be reconstituted as follow :  
- mab0050-0 in 0,5 ml of distilled water  
- mab0050-1 in 1 ml of distilled water

**Preservative** None

**Storage** Short term storage : +4°C  
Long term storage : -20°C

#### References

Muzerelle A et al; Tetanus neurotoxin-insensitive vesicle-associated membrane protein localizes to a presynaptic membrane compartment in selected terminal subsets of the rat brain. Neuroscience. 2003;122:59-75

#### Product type

Primary antibody

#### Clonality

Monoclonal Antibody

#### Related product

pab0053 ; pab0054 ;  
pab0056 ; pab0035 ;  
lab0235 ; lab0239 ;  
lab0242 ; lab0246 ;  
lab0250 ; lab0252 ;  
lab0257 ; opr0002 ;  
opr0008 ; opr0009

#### Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated

version A

**Covalab France**  
11, Avenue Albert Einstein  
69100 Villeurbanne - France  
Phone +33 (0) 437 654 230  
Fax +33 (0) 437 289 416  
Email [pa@covalab.com](mailto:pa@covalab.com)  
[www.covalab.com](http://www.covalab.com)

**Covalab - UK Ltd**  
St John's Innovation Centre  
Unit 75, Cowley Road  
Cambridge, CB40WS - UK  
Phone +44 (0) 1223 421055  
Fax +44 (0) 1223 420844  
Email [enquiries@covalab.co.uk](mailto:enquiries@covalab.co.uk)  
[www.covalab.co.uk](http://www.covalab.co.uk)