

**Opr0006 Covabtest® NHS Activated Plates**

**Principle**

These functionalized plates allow macromolecules or haptens containing free primary amines to be covalently immobilised.

**Application**

Determining antibody levels in the serum of immunized animal or in cell supernatant from hybridoma clones, Affinity measurements by competition Elisa, Adhesion-molecule assays, Enzyme-substrate studies, Ligand-receptor interactions...

**Packaging instructions**

**Packaging** 5 Plates  
**Appearance** Plates, Blocking solution 100 mL  
**Storage instructions** Store at +4°C

**Product type**

Plates

**Related product(s)**

opr0004 ; opr0005

Updated 19/02/2004  
Version A

**Covalab France**  
Bio-Pôle Lyon Gerland  
3, Allée Pierre de Coubertin  
69007 Lyon - France  
Phone +33 (0) 437 654 230  
Fax +33 (0) 437 289 416  
Email [pa@covalab.com](mailto:pa@covalab.com)  
[www.covalab.com](http://www.covalab.com)

**Covalab - UK Ltd**  
St John's Innovation Centre  
Unit 75, Cowley Road  
Cambridge, CB40WS - UK  
Phone +44 (0) 1223 421055  
Fax +44 (0) 1223 420844  
Email [enquiries@covalab.co.uk](mailto:enquiries@covalab.co.uk)  
[www.covalab.co.uk](http://www.covalab.co.uk)