

Opr0011 Streptavidin Peroxydase Conjugate

Product specifications

Working dilution

Working dilution should be determined by titration assay.



The following are guidelines only

- ELISA: 1/2000 to 1/20 000
- Western Blot: 1/500 to 1/2000
- Immunohistology: 1/500 to 1/5000

Reactivity

By the ELISA method using polystyrene wells coated with biotinylated BSA

By chemiluminescence after PAGE & transfer to nitrocellulose of biotinylated molecular weigh markers (Ref opr0002)



Product type

Other product

Related product(s)

Packaging specifications

Packaging 0.5 ml

Appearance Liquid (Clear colourless solution)

Constituents PBS pH 7.4, 1% BSA, 50% Glycerol

Preservative 0.02% Thimerosal

Storage instructions For continuous use, store at 2-8°C for up to one month
For extended storage, the solution may be stored at -20°C

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

Updated Version 07/06/2005 A

Covalab France
Bio-Pôle Lyon Gerland
3, Allée Pierre de Coubertin
69007 Lyon – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd
St John’s Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk