

## PRODUCT DATA SHEET

### pab0007 Polyclonal Antibody to Maf protein



#### Product Number/Size

pab0007-0 Lyophilized from 100 µl of serum  
pab0007-1 Lyophilized from 200 µl of serum

#### Product specifications

**Immunogen** Maf-GST fusion protein

**Form** Antiserum **Raised in** Rabbit

**Specificity** Antiserum to avian Maf protein reacts specifically with 47 kDa protein derived from a transfected CEFcell line

**Cross reactions** Not yet tested

**Working dilutions** Optimal dilution should be determined by the end user  
The following are guidelines only :  
- Immunohistology : 1/50 to 1/500  
- Western Blot : 1/200 to 1/5 000  
Useful for electrophoretic mobility shift assay (EMSA)

#### Packaging specifications

**Appearance** Lyophilized powder (glass vial)  
Vial must be reconstituted as follow :  
- pab0007-0 in 100 µl of distilled water  
- pab0007-1 in 200 µl of distilled water

**Preservative** None

**Storage** Short term storage : +4°C  
Long term storage : -20°C

#### References

#### Product type

Primary antibody

#### Clonality

Polyclonal Antibody

#### Related product

lab0273 ; lab0270 ;  
lab0264 ; lab0271 ;  
lab0274 ; lab0269 ;  
lab0272 ; opr0002 ;  
opr0008 ; opr0009

#### Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 08/02/2008

version C

#### Covalab France

11, Avenue Albert Einstein  
69100 Villeurbanne – France  
Phone +33 (0) 437 654 230  
Fax +33 (0) 437 289 416  
Email [pa@covalab.com](mailto:pa@covalab.com)  
[www.covalab.com](http://www.covalab.com)

#### Covalab – UK Ltd

St John's Innovation Centre  
Unit 75, Cowley Road  
Cambridge, CB40WS – UK  
Phone +44 (0) 1223 421055  
Fax +44 (0) 1223 420844  
Email [enquiries@covalab.co.uk](mailto:enquiries@covalab.co.uk)  
[www.covalab.co.uk](http://www.covalab.co.uk)