

PRODUCT DATA SHEET

pab0025 Polyclonal Antibody to Phosphothreonine

Product Number/Size

pab0025-0 Lyophilized from 100 µl of serum
pab0025-1 Lyophilized from 200 µl of serum

Product specifications

Immunogen	KLH-Phosphothreonine		
Form	Antiserum	Raised in	Rabbit
Specificity	Antiserum reacts with phosphorylated threonine both as free amino acid or when conjugated to a carrier protein such as BSA		
Cross reactions	Cross reacts with phosphorylated tyrosine and does not cross react with phosphorylated serine		
Working dilutions	Optimal dilution should be determined by the end user The following are guidelines only : - Immunohistology : 1/50 to 1/500 - Western Blot : 1/200 to 1/2 000		

Packaging specifications

Appearance	Lyophilized powder (glass vial) Vial must be reconstituted as follow : - pab0025-0 in 100 µl of distilled water - pab0025-1 in 200 µl of distilled water
Preservative	None
Storage	Short term storage : +4°C Long term storage : -20°C

References



Product type

Primary antibody

Clonality

Polyclonal Antibody

Related product

lab0273 ; lab0270 ;
lab0264 ; lab0271 ;
lab0274 ; lab0269 ;
lab0272 ; opr0002 ;
opr0008 ; opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 08/02/2008

version C

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk