

## PRODUCT DATA SHEET

### pab0031 Polyclonal Antibody to Human Erythrocyte Band EP4.2 protein

#### Product Number/Size

pab0031-0 Lyophilized from 100 µl of serum  
pab0031-1 Lyophilized from 200 µl of serum

#### Product specifications

<b>Immunogen</b>	A 15 residue synthetic peptide derived from specific domain of band EP4.2 protein		
<b>Form</b>	Antiserum	<b>Raised in</b>	Rabbit
<b>Specificity</b>	Antiserum reacts specifically with band EP4.2 protein from human erythrocyte cells		
<b>Cross reactions</b>	Does not cross react with transglutaminases		
<b>Working dilutions</b>	Optimal dilution should be determined by the end user The following are guidelines only : - Immunohistochemistry : 1/200 to 1/1 000 (For IHC Immunopurified Ab is recommended: see related products) - Western Blot : 1/200 to 1/1 000		

#### Packaging specifications

<b>Appearance</b>	Lyophilized powder (glass vial) Vial must be reconstituted as follow : - pab0031-0 in 100 µl of distilled water - pab0031-1 in 200 µl of distilled water
<b>Preservative</b>	None
<b>Storage</b>	Short term storage : +4°C Long term storage : -20°C

#### References



#### Product type

Primary antibody

#### Clonality

Polyclonal Antibody

#### Related product

pab0032 ; lab0273 ;  
lab0270 ; lab0264 ;  
lab0271 ; lab0274 ;  
lab0269 ; lab0272 ;  
opr0002 ; opr0008 ;  
opr0009

#### Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 08/02/2008

version C

#### Covalab France

11, Avenue Albert Einstein  
69100 Villeurbanne – France  
Phone +33 (0) 437 654 230  
Fax +33 (0) 437 289 416  
Email [pa@covalab.com](mailto:pa@covalab.com)  
[www.covalab.com](http://www.covalab.com)

#### Covalab – UK Ltd

St John's Innovation Centre  
Unit 75, Cowley Road  
Cambridge, CB40WS – UK  
Phone +44 (0) 1223 421055  
Fax +44 (0) 1223 420844  
Email [enquiries@covalab.co.uk](mailto:enquiries@covalab.co.uk)  
[www.covalab.co.uk](http://www.covalab.co.uk)