

PRODUCT DATA SHEET

pab0034 Immunopurified Antibody to Prion Protein PrP

Product Number/Size

pab0034 25 µg of purified antibody

Product specifications

Immunogen	A 15 residue synthetic peptide derived from C-terminal domain of bovin PrP protein		
Form	Peptide affinity purified	Raised in	Rabbit
Specificity	Purified antibody reacts specifically with full sequence PrP protein of brain tissue from different species (bovin, human, sheep, rat ...)		

Cross reactions

Working dilutions Optimal dilution should be determined by the end user
The following are guidelines only :
- Immunohistochemistry : 1/20 to 1/100
- Western Blot : 1/20 to 1/100

Packaging specifications

Appearance	Liquid (glass vial)
Preservative	0.02% NaN ₃
Storage	Stored at -20°C.

References



Product type

Primary antibody

Clonality

Polyclonal Antibody

Related product

pab0033 ; mab0026 ;
lab0273 ; lab0270 ;
lab0264 ; lab0271 ;
lab0274 ; lab0269 ;
lab0272 ; opr0002 ;
opr0008 ; opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 04/10/06

version B

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk