

PRODUCT DATA SHEET

pab0040 Polyclonal Antibody to Human Caspase-3



Product Number/Size

pab0040-0 Lyophilized from 100 µl of serum
pab0040-1 Lyophilized from 200 µl of serum

Product specifications

Immunogen	Epitope corresponding to an amino acid sequence mapping from C-terminal domain of human cysteine protease CPP32 (Caspase-3)	
Form	Antiserum	Raised in Rabbit
Specificity	Antiserum to Caspase-3 (CPP32) reacts with 32 kDa protein (Caspase-3) derived from a HL60 and U937 cells	
Cross reactions	Not yet tested	
Working dilutions	Optimal dilution should be determined by the end user The following are guidelines only : - Immunohistology : 1/100 to 1/500 - Western Blot : 1/500 to 1/5 000	

Packaging specifications

Appearance	Lyophilized powder (glass vial) Vial must be reconstituted as follow : - pab0040-0 in 100 µl of distilled water - pab0040-1 in 200 µl of distilled water
Preservative	None
Storage	Short term storage : +4°C Long term storage : -20°C

References

Product type

Primary antibody

Clonality

Polyclonal Antibody

Related product

pep0010 ; lab0273 ;
lab0270 ; lab0264 ;
lab0271 ; lab0274 ;
lab0269 ; lab0272 ;
opr0002 ; opr0008 ;
opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 08/02/2008

version C

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk