

## PRODUCT DATA SHEET

### pab0044 Immunopurified Antibody to Peptidylarginine deiminase type 1 (PAD1)

#### Product Number/Size

pab0044 25 µg of purified antibody

#### Product specifications

<b>Immunogen</b>	A 15 residue synthetic peptide derived from a domaine specific for human peptidylarginine deiminase type 1 protein		
<b>Form</b>	Peptide affinity purified	<b>Raised in</b>	Rabbit
<b>Specificity</b>	Purified antibody reacts with peptidylarginine deiminase type 1 (PAD1)		
<b>Cross reactions</b>	Reactivity with other PAD isoforms has not yet been determined		
<b>Working dilutions</b>	Optimal dilution should be determined by the end user The following are guidelines only : - Immunohistochemistry : 1/20 to 1/100 - Western Blot : 1/20 to 1/200		

#### Packaging specifications

<b>Appearance</b>	Liquid (glass vial)
<b>Preservative</b>	0.02% NaN <sub>3</sub>
<b>Storage</b>	Stored at -20°C.

#### References



#### Product type

Primary antibody

#### Clonality

Polyclonal Antibody

#### Related product

pab0043 ; lab0273 ;  
lab0270 ; lab0264 ;  
lab0271 ; lab0274 ;  
lab0269 ; lab0272 ;  
opr0002 ; opr0008 ;  
opr0009

#### Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 04/10/06

version B

#### Covalab France

11, Avenue Albert Einstein  
69100 Villeurbanne – France  
Phone +33 (0) 437 654 230  
Fax +33 (0) 437 289 416  
Email [pa@covalab.com](mailto:pa@covalab.com)  
[www.covalab.com](http://www.covalab.com)

#### Covalab – UK Ltd

St John's Innovation Centre  
Unit 75, Cowley Road  
Cambridge, CB40WS – UK  
Phone +44 (0) 1223 421055  
Fax +44 (0) 1223 420844  
Email [enquiries@covalab.co.uk](mailto:enquiries@covalab.co.uk)  
[www.covalab.co.uk](http://www.covalab.co.uk)