

PRODUCT DATA SHEET

pab0061 Immunopurified Antibody to Transglutaminase 1

Product Number/Size

pab0061 25 µg of purified antibody

Product specifications

Immunogen	A 15 residue synthetic peptide derived from specific domain of TGase 1 protein		
Form	Peptide affinity purified	Raised in	Rabbit
Specificity	Purified antibody reacts specifically with keratinocyte transglutaminase 1 from human		
Cross reactions	Not yet tested		
Working dilutions	Optimal dilution should be determined by the end user The following are guidelines only : - Immunohistochemistry : 1/20 to 1/100 - Western Blot : 1/20 to 1/200		

Packaging specifications

Appearance	Liquid (glass vial)
Preservative	0.02% NaN ₃
Storage	Stored at -20°C.

References



Product type

Primary antibody

Clonality

Polyclonal Antibody

Related product

pab0060 ; lab0273 ;
lab0270 ; lab0264 ;
lab0271 ; lab0274 ;
lab0269 ; lab0272 ;
opr0002 ; opr0008 ;
opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 04/10/06

version B

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk