

## PRODUCT DATA SHEET

### pab0062 Polyclonal Antibody to Transglutaminase 2

#### Product Number/Size

pab0062-0 Lyophilized from 100 µl of serum  
pab0062-1 Lyophilized from 200 µl of serum

#### Product specifications

**Immunogen** A 15 residue synthetic peptide derived from specific domain of human TGase 2 protein (tTG)

**Form** Antiserum **Raised in** Rabbit

**Specificity** Antiserum to TGase 2 reacts specifically with TGase 2 protein from human cells

**Cross reactions** No cross reactivity with guinea pig TGase 2 neither with any TGase isoformes

**Working dilutions** Optimal dilution should be determined by the end user  
The following are guidelines only :  
- Immunohistochemistry : 1/200 to 1/2 000 (For IHC Immunopurified Ab is recommended: see related products)  
- Western Blot : 1/200 to 1/1 000

#### Packaging specifications

**Appearance** Lyophilized powder (glass vial)  
Vial must be reconstituted as follow :  
- pab0062-0 in 100 µl of distilled water  
- pab0062-1 in 200 µl of distilled water

**Preservative** None

**Storage** Short term storage : +4°C  
Long term storage : -20°C

#### References



#### Product type

Primary antibody

#### Clonality

Polyclonal Antibody

#### Related product

pab0063 ; lab0273 ;  
lab0270 ; lab0264 ;  
lab0271 ; lab0274 ;  
lab0269 ; lab0272 ;  
opr0002 ; opr0008 ;  
opr0009

#### Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 05/02/08

version C

#### Covalab France

11, Avenue Albert Einstein  
69100 Villeurbanne – France  
Phone +33 (0) 437 654 230  
Fax +33 (0) 437 289 416  
Email [pa@covalab.com](mailto:pa@covalab.com)  
[www.covalab.com](http://www.covalab.com)

#### Covalab – UK Ltd

St John's Innovation Centre  
Unit 75, Cowley Road  
Cambridge, CB40WS – UK  
Phone +44 (0) 1223 421055  
Fax +44 (0) 1223 420844  
Email [enquiries@covalab.co.uk](mailto:enquiries@covalab.co.uk)  
[www.covalab.co.uk](http://www.covalab.co.uk)