

PRODUCT DATA SHEET

pab0066 Polyclonal Antibody to Transglutaminase 4

Product Number/Size

pab0066-0 Lyophilized from 100 µl of serum
pab0066-1 Lyophilized from 200 µl of serum

Product specifications

Immunogen	A 15 residue synthetic peptide derived from specific domain of Tgase 4 protein		
Form	Antiserum	Raised in	Rabbit
Specificity	Antiserum to TGase 4 reacts specifically with prostatic transglutaminase 4 (TGase 4) from human prostate cell		
Cross reactions	Not yet tested		
Working dilutions	Optimal dilution should be determined by the end user The following are guidelines only : - Immunohistochemistry : 1/200 to 1/500 (For IHC Immunopurified Ab is recommended: see related products) - Western Blot : 1/200 to 1/1 000		

Packaging specifications

Appearance	Lyophilized powder (glass vial) Vial must be reconstituted as follow : - pab0066-0 in 100 µl of distilled water - pab0066-1 in 200 µl of distilled water
Preservative	None
Storage	Short term storage : +4°C Long term storage : -20°C

References



Product type

Primary antibody

Clonality

Polyclonal Antibody

Related product

pab0067 ; lab0273 ;
lab0270 ; lab0264 ;
lab0271 ; lab0274 ;
lab0269 ; lab0272 ;
opr0002 ; opr0008 ;
opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 05/02/08

version C

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk