

PRODUCT DATA SHEET

pab0218 Polyclonal Antibody to Human Podoplanin

Product Number/Size

pab0218-0 Lyophilized from 100 µl of serum
pab0218-1 Lyophilized from 200 µl of serum

Product specifications

Immunogen Synthetic peptide derived from C terminal region of human Podoplanin

Form Antiserum **Raised in** Rabbit

Specificity Antiserum reacts specifically with 25 kDa protein

Cross reactions Not yet tested

Working dilutions Optimal dilutions should be determined by the end user
The following are guidelines only :
- FACS : 1/50 to 1/200
- Immunohistochemistry : 1/50 to 1/500
- Western Blot : 1/200 to 1/1 000

Packaging specifications

Appearance Lyophilized powder (glass vial)
Vial must be reconstituted as follow :
- pab0218-0 in 100 µl of distilled water
- pab0218-1 in 200 µl of distilled water

Preservative None

Storage Short term storage : +4°C
Long term storage : -20°C

References

Nisato R et al; Generation and characterization of telomerase-transfected human lymphatic endothelial cells with an extended life span. Am J Pathol. 2004; 165(1):11-24.



Product type

Primary antibody

Clonality

Polyclonal antibody

Related product

lab0270; lab0264;
lab0271; lab0274;
lab0269; lab0272;
lab0273; opr0002;
opr0008; opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 08/02/2008

version B

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne - France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab - UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS - UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk