

PRODUCT DATA SHEET

pab0226 Polyclonal Antibody to E3 ubiquitin-protein ligase (SIAH-2)

Product Number/Size

pab0226-0 Lyophilized from 100 µl of serum
pab0226-1 Lyophilized from 200 µl of serum

Product specifications

Immunogen Synthetic peptide derived from N terminal region of SIAH-2

Form Antiserum **Raised in** Rabbit

Specificity Antiserum reacts specifically with human 35 kDa SIAH-2 protein

Cross reactions Cross reacts with all species

Working dilutions Optimal dilution should be determined by the end user
The following are guidelines only :
- Immunohistochemistry : 1/100 to 1/500
- Western Blot : 1/500 to 1/2 000

Packaging specifications

Appearance Lyophilized powder (glass vial)
Vial must be reconstituted as follow :
- pab0226-0 in 100 µl of distilled water
- pab0226-1 in 200 µl of distilled water

Preservative None

Storage Short term storage : +4°C
Long term storage : -20°C

References

Fiucci G. et al; Siah-1b is a direct transcriptional target of p53: identification of the functional p53 responsive element in the siah-1b promoter.
Proc Natl Acad Sci U S A. 2004;101(10):3510-5.



Product type

Primary antibody

Clonality

Polyclonal antibody

Related product

lab0270; lab0264;
lab0271; lab0274;
lab0269; lab0272;
lab0273; opr0002;
opr0008; opr0009

Warning

This is a laboratory reagent. It is not to be administered to human or animals nor be used as a drug

updated 05/02/08

version B

Covalab France

11, Avenue Albert Einstein
69100 Villeurbanne – France
Phone +33 (0) 437 654 230
Fax +33 (0) 437 289 416
Email pa@covalab.com
www.covalab.com

Covalab – UK Ltd

St John's Innovation Centre
Unit 75, Cowley Road
Cambridge, CB40WS – UK
Phone +44 (0) 1223 421055
Fax +44 (0) 1223 420844
Email enquiries@covalab.co.uk
www.covalab.co.uk